

Internal Hexagon Prosthetic Options

MIS[®] | MAKE IT SIMPLE

NOT FOR USE IN THE USA AND CANADA

Products may not be available in some regulated markets. Please contact your local Dentsply Sirona sales office for current product assortment and availability.

The MIS Quality System complies with international quality standard ISO 13485: 2016 - Quality Management System for Medical Devices, with Medical Device Directive 93/42/EEC, and with EU Medical Device Regulation MDR 2017/745. Please note, that not all products are registered or available in every country/region.

IFUs for MIS products may be found at: <https://ifu.mis-implants.com>. Adobe Acrobat is required to view the IFU file on the website. This software may be freely downloaded from the Adobe website.

Table of Contents

Implant Dimensions P. 6 - 7	Prosthetic Options Listed by Healing Caps P. 8 - 9	Healing Caps P. 10 - 12	Impression Coping P. 13 - 14	Temporary Restoration P. 15 - 17
---------------------------------------	--	-----------------------------------	--	--

○ Temporary Cylinders

Screw-Retained Restoration

P. 18 - 23

- Multi-Unit System
- Plastic Cylinders

Digital Restoration

P. 24 - 29

- EZ-Base System
- Ti-Bases
- Digital Model Analogs
- Scan Posts
- Titanium Blanks

Cemented Restoration

P. 30 - 39

- Cementable Abutments
- CPK

Overdenture

P. 40 - 45

- LOCKIT System
- OT-Equator Kits
- Ball Attachment System

Prosthetic Screws and Drivers

P. 46 - 47

Implant Dimensions: SEVEN





Narrow
Platform



Standard
Platform



Wide
Platform





Implant Dimensions: M4




Narrow
Platform



Standard
Platform



Wide
Platform




ALL PLATFORMS

Prosthetic Options Listed by Healing Caps


Each column lists the different heights of healing caps and corresponding prosthetic options. **Bold** item numbers indicate the products shown in the images above.

*Image is for illustration only. (2mm height closed tray impression coping presented. Not available in 3mm).


1. Narrow Platform



2. Standard Platform




3. Wide Platform



Healing caps	Impression coping, open tray	Impression coping, closed tray	Cemented options	Digital options
--------------	------------------------------	--------------------------------	------------------	-----------------

or

or




Gingival heights 2, 3, 4mm	MH-N0440	MN-IO401	MN-ITC30*	MK-NPC62 MN-MC102 MN-AT102 MN-AT202	
	MH-N0540			MK-NPC63	MN-TBC30 MN-FTB30
	MH-N0640			MK-NPC64	

Healing caps	Impression coping, open tray	Impression coping, closed tray	Screw retained options	Cemented options	Digital options
--------------	------------------------------	--------------------------------	------------------------	------------------	-----------------


or or or

Healing caps	Impression coping, open tray	Impression coping, closed tray	Screw retained options	Cemented options	Digital options
--------------	------------------------------	--------------------------------	------------------------	------------------	-----------------

or or or



Gingival heights 2, 3, 4mm	MH-S0448	MN-IO481	MD-ITC30*	MK-CPC62 MD-MC102 MD-AT152 MD-AT252	MA-SC248	
	MH-S0548			MD-AT153 MD-AT253 MK-CPC63	MA-SC348	MD-TBC30 MD-FTB30
	MH-S0648			MK-CPC64	MA-SC448	



Gingival heights 2, 3, 4mm	MH-W0455	MN-IO551	MW-ITC30*	MK-WPC62 MW-MC102 MW-AT152		
	MH-W0555			MK-WPC63	MA-WC348	MW-TBC30 MW-FTB30
	MH-W0655			MK-WPC64		

Healing Caps Narrow Platform

Tightening torque: 20-25 Ncm

HEIGHT

SEVEN/M4
H2

SEVEN/M4
H3

SEVEN/M4
H4

SEVEN/M4
H5

SEVEN/M4
H6

SEVEN/M4
H8

Thin, Ø3.3



MH-N2330




MH-N3330



MH-N4330




MH-N6330




MH-N8330


Standard, Ø4



MH-N0440



MH-N0540



MH-N0640


Anatomic, Ø4.3



MH-52330




MH-53330



MH-54330



MH-55330



MH-56330

Healing Caps Standard Platform

Tightening torque: 20-25 Ncm

HEIGHT

SEVEN/M4
H3

SEVEN/M4
H4

SEVEN/M4
H5

SEVEN/M4
H6

Thin, Ø4



MH-03375



MH-04375



MH-05375



MH-06375

Standard, Ø4.8



MH-A3375



MH-A4375



MH-A5375



MH-A6375

Standard, Ø4.8



MH-S0448



MH-S0548




MH-S0648


Anatomic, Ø5.5




MH-53375



MH-54375



MH-55375




MH-56375

Healing Caps Narrow Platform


Tightening torque: 20-25 Ncm

HEIGHT


Standard, Ø5



Standard, Ø5.5



Anatomic, Ø6.5



SEVEN/M4
H3

SEVEN/M4
H4

SEVEN/M4
H5

SEVEN/M4
H6

Impression Copings Open tray

Tightening torque: 20-25 Ncm




PLATFORM

SEVEN/M4
NP




SEVEN/M4
SP

SEVEN/M4
WP




Impression coping



Impression coping



Implant analog,
open/closed tray



Impression Copings

Closed tray

Tightening torque: 20-25 Ncm




PLATFORM

SEVEN/M4
NP






SEVEN/M4
SP

SEVEN/M4
WP



Impression coping



Impression coping set



Press fit



Temporary Restoration

Temporary Plastic Cylinders

Anti rotation

Tightening torque: 20-25 Ncm


PLATFORM

Narrow profile temporary plastic abutment, H1

Temporary plastic abutment, H1


Temporary plastic cylinder, H2.5

SEVEN/M4
NP




MN-TPA12

SEVEN/M4
SP




MD-TPA12


SEVEN/M4
WP




MW-TPA12




MN-TPA10




MD-TPA10



MW-TPA10



MD-TPH50



MW-TPH50

Temporary Restoration

Temporary Cylinders

Free rotation

Tightening torque: 20-25 Ncm



PLATFORM

SEVEN/M4
NP



SEVEN/M4
SP

SEVEN/M4
WP



H1



H2



H3



Temporary Restoration

Temporary Cylinders

Anti-rotation

Tightening torque: 20-25 Ncm




PLATFORM

SEVEN/M4
NP



SEVEN/M4
SP

SEVEN/M4
WP




H1



H2



H3



Screw Retained Restoration

Multi-Unit System

Tightening torque: 30 Ncm

PLATFORM

SEVEN/M4
NP

SEVEN/M4
SP

SEVEN/M4
WP

H2



MA-NC248



MA-SC248

H3



MA-NC348



MA-SC348



MA-WC348

H4



MA-NC448



MA-SC448

H5



MA-WC548

Screw Retained Restoration

Multi-Unit System

Tightening torque: 30 Ncm

PLATFORM

SEVEN/M4
NP

SEVEN/M4
SP

SEVEN/M4
WP

H1



MA-N1480



MA-S1480



MA-W1480

H2



MA-N2480



MA-S2480

H3



MA-N3480



MA-S3480




MA-W3480

H4




MA-N4480



MA-S4480


H5



MA-N5480



MA-S5480



MA-W5480

Screw Retained Restoration

Multi-Unit System

Angulated abutments, 17°

Tightening torque: 20 Ncm

PLATFORM

SEVEN/M4
NP

SEVEN/M4
SP

SEVEN/M4
WP

H1



MA-NO171



MA-SO171



MA-WO171

H2



MA-NO172



MA-SO172



MA-WO172

H3



MA-SO173



MA-WO173

Screw Retained Restoration

Multi-Unit System

Angulated abutments, 30°

Tightening torque: 20 Ncm

PLATFORM

SEVEN/M4
NP

SEVEN/M4
SP

SEVEN/M4
WP

H1



MA-NO301




MA-SO301




MA-WO301


H2



MA-NO302




MA-SO302




MA-WO302

H3




MA-SO303



MA-WO303


Multi-Unit System Components

Healing cap




MU-H4480

Analog



MU-RSM48

Impression coping,
open tray



MU-I0480

Impression coping,
closed tray



MU-IT480

Temporary cylinder



MU-TO480

Plastic cylinder



MU-CP048

Gold bar coping abutment



MU-BGC48

Cementing cap
(New design)



MU-CC481

Cementing cap



MU-CC480

Scan Post



MU-SP102

Digital model analog



MU-MAN10

Screw Retained Restoration

Plastic Cylinders

Free rotation

Tightening torque: 30 Ncm


PLATFORM

SEVEN/M4
NP



SEVEN/M4
SP

SEVEN/M4
WP



Gold plastic cylinder





Plastic cylinder



CoCr plastic cylinder



Wide profile plastic cylinder



Screw Retained Restoration
Plastic Cylinders
 Anti-rotation

Tightening torque: 30 Ncm




PLATFORM

SEVEN/M4
 NP




SEVEN/M4
 SP

SEVEN/M4
 WP




Gold plastic cylinder





Plastic cylinder



CoCr plastic cylinder



Wide profile plastic cylinder



Digital Restoration

EZ-Base System

Tightening torque: 25 Ncm

PLATFORM

H1.5



SEVEN/M4
NP



SEVEN/M4
SP




SEVEN/M4
WP




EZ-Base System Components

Ratchet key for angled screw channel




MT-ELR10 - Long
MT-ESR10 - Short

Motor key for angled screw channel



MT-ELM10 - Long
MT-ESM10 - Short

EZ-Base prosthetic screw



MN-E0160 - NP
MD-E0220 - SP



A unique screw and driver tip design enable extreme angulation and a wide insertion path.

Digital Restoration

Ti-Base

Anti-rotation

Tightening torque: 30 Ncm



PLATFORM

SEVEN/M4
NP




SEVEN/M4
SP

SEVEN/M4
WP



H0.5





H1.5



H3



H0.5



Digital Restoration

Ti-Base

Free-rotation

Tightening torque: 30 Ncm



PLATFORM

SEVEN/M4
NP




SEVEN/M4
SP

SEVEN/M4
WP



H0.5



H1.5



H3



Digital Restoration
**Digital Model
Analog**

Tightening torque: 15 Ncm


SEVEN/M4
NP




SEVEN/M4
SP



SEVEN/M4
WP



Multi-Unit



Digital Restoration
Scan Posts

Tightening torque: 20 Ncm

SEVEN/M4
NP



MN-SP102

SEVEN/M4
SP



MD-SP102

SEVEN/M4
WP




MW-SP102

Multi-Unit



MU-SP102



MK-0057
SCAN POST KIT,
INTERNAL HEX.

Digital Restoration
Titanium Blanks
 Anti-rotation

Tightening torque: 30 Ncm

PLATFORM

SEVEN/M4
 NP

SEVEN/M4
 SP

SEVEN/M4
 WP

Ø9

Ø9mm



MN-IB091

Ø9mm



MD-IB091

Ø9mm



MW-IB091

Ø12

Ø12mm



MN-IB121

Ø12mm



MD-IB121

Ø12mm



MW-IB121

Amann Girrbach, Ø12

Ø12mm



MN-AGT12

Ø12mm



MD-AGT12

Ø12mm



MW-AGT12

Milling analogs

Ø5



MN-MAB10

13.5mm

Ø5



MD-MAB10

13.5mm

Ø5



MW-MAB10

13.5mm



Cemented Restoration
**Cementable
Abutments**

Tightening torque: 30 Ncm



PLATFORM

SEVEN/M4
NP

10°



20°




Cemented Restoration
**Cementable
 Abutments**

Tightening torque: 30 Ncm

SEVEN/M4
 SP

PLATFORM



Cemented Restoration


Cementable Abutments

Tightening torque: 30 Ncm




PLATFORM

SEVEN/M4
WP

0°



15°



Cemented Restoration
Cementing posts

Tightening torque: 30 Ncm


PLATFORM

SEVEN/M4
NP


SEVEN/M4
SP

SEVEN/M4
WP


H1



H1



H2



Cemented Restoration
**CPK -
 Transgingival
 Abutments**
 4mm Crown

Tightening torque: 30 Ncm




PLATFORM

SEVEN/M4
NP




SEVEN/M4
SP

SEVEN/M4
WP




H1






H2



H3



H4



Cemented Restoration
**CPK -
 Transgingival
 Abutments**
 6mm Crown

Tightening torque: 30 Ncm


PLATFORM

SEVEN/M4
 NP


SEVEN/M4
 SP

SEVEN/M4
 WP


H1



MN-CPK61




MD-CPK61




MW-CPK61


H2



MN-CPK62




MN-CPC62




MD-CPK62



MD-CPC62




MW-CPK62



MW-CPC62


H3



MN-CPK63




MN-CPC63




MD-CPK63



MD-CPC63




MW-CPK63



MW-CPC63


H4




MN-CPK64




MN-CPC64




MD-CPK64



MD-CPC64



MW-CPK64



MW-CPC64

Cemented Restoration

CPK - Transgingival Abutments

8mm Crown


Tightening torque: 30 Ncm

PLATFORM


SEVEN/M4
SP

SEVEN/M4
WP

H1




MD-CPK81




MW-CPK81

H2




MD-CPK82




MW-CPK82

H3




MD-CPK83




MW-CPK83

H4



MD-CPK84



MW-CPK84

CPK abutments are also sold in kits.

MK-NPK41 ● Ø4, H=1, C=4
MK-CPK41 ● Ø4.8, H=1, C=4
MK-WPK41 ● Ø5.5, H=1, C=4

MK-NPK42 ● Ø4, H=2, C=4
MK-CPK42 ● Ø4.8, H=2, C=4
MK-WPK42 ● Ø5.5, H=2, C=4

MK-NPK43 ● Ø4, H=3, C=4
MK-CPK43 ● Ø4.8, H=3, C=4
MK-WPK43 ● Ø5.5, H=3, C=4

MK-NPK44 ● Ø4, H=4, C=4
MK-CPK44 ● Ø4.8, H=4, C=4
MK-WPK44 ● Ø5.5, H=4, C=4

MK-NPK61 ● Ø4, H=1, C=6
MK-CPK61 ● Ø4.8, H=1, C=6
MK-WPK61 ● Ø5.5, H=1, C=6

MK-NPC62 ● Ø4, H=2, C=6
MK-CPC62 ● Ø4.8, H=2, C=6
MK-WPC62 ● Ø5.5, H=2, C=6

MK-NPK62 ● Ø4, H=2, C=6
MK-CPK62 ● Ø4.8, H=2, C=6
MK-WPK62 ● Ø5.5, H=2, C=6

MK-NPC63 ● Ø4, H=3, C=6
MK-CPC63 ● Ø4.8, H=3, C=6
MK-WPC63 ● Ø5.5, H=3, C=6

MK-NPK63 ● Ø4, H=3, C=6
MK-CPK63 ● Ø4.8, H=3, C=6
MK-WPK63 ● Ø5.5, H=3, C=6

MK-NPC64 ● Ø4, H=4, C=6
MK-CPC64 ● Ø4.8, H=4, C=6
MK-WPC64 ● Ø5.5, H=4, C=6

MK-NPK64 ● Ø4, H=4, C=6
MK-CPK64 ● Ø4.8, H=4, C=6
MK-WPK64 ● Ø5.5, H=4, C=6

MK-CPK81 ● Ø4.8, H=1, C=8
MK-WPK81 ● Ø5.5, H=1, C=8

MK-CPK82 ● Ø4.8, H=2, C=8
MK-WPK82 ● Ø5.5, H=2, C=8

MK-CPK83 ● Ø4.8, H=3, C=8
MK-WPK83 ● Ø5.5, H=3, C=8

MK-CPK84 ● Ø4.8, H=4, C=8
MK-WPK84 ● Ø5.5, H=4, C=8



Burnout plastic cap,
anti-rotation




MN-ICH40
MD-ICH40
MW-ICH40

Burnout plastic cap,
free rotation




MN-IC040
MD-IC040
MW-IC040

Plastic healing cap




MM-CHC40
MM-CHC60
MM-CHC80

Transfer plastic cap



MN-IC600
MD-IC800
MW-IC800

Analogs



MN-RSM60
MD-RSM60
MW-RSM60

Overdenture
LOCKiT
 Attachment
 System

Tightening torque: 30 Ncm

i All LOCKiT products shown in this catalog are sold subject to individual regulatory certification only.





PLATFORM

SEVEN/M4
 NP

SEVEN/M4
 SP

SEVEN/M4
 WP

H1




MW-LOC01

H2

MN-LOC02




MD-LOC02



H3

MN-LOC03



MD-LOC03



MW-LOC03

H4

MN-LOC04



MD-LOC04



H5

MN-LOC05

MD-LOC05

MW-LOC05

LOCKiT System Components



LOCKiT abutments
are also sold in kits

- MK-NLC1 ● Ø3.75, H=1
- MK-SLC1 ● Ø3.75, H=1
- MK-WLC1 ● Ø4.5, H=1
- MK-NLC2 ● Ø3.75, H=2
- MK-SLC2 ● Ø3.75, H=2
- MK-NLC3 ● Ø3.75, H=3
- MK-SLC3 ● Ø3.75, H=3
- MK-WLC3 ● Ø4.5, H=3
- MK-NLC4 ● Ø3.75, H=4
- MK-SLC4 ● Ø3.75, H=4
- MK-NLC5 ● Ø3.75, H=5
- MK-SLC5 ● Ø3.75, H=5
- MK-WLC5 ● Ø4.5, H=5

Core tool



MT-LOC10

Ratchet driver



MT-RSD19 - Short
MT-RSD24 - Long

Caps range



MM-BLA10
MM-RED10
MM-GRE10
MM-BLU10
MM-PIN10
MM-CLE10

Male processing
package



MM-LOC02

Metal housing



MM-TIT10

Impression coping



MM-LOC11

Analog



MM-LOA40 - NP/SP
MM-LOA50 - WP

Abutment holder sleeve



MT-LOC11

Overdenture
OT-Equators

Tightening torque: 30 Ncm

PLATFORM

SEVEN/M4
NP

SEVEN/M4
SP

SEVEN/M4
WP

H1



MK-NOE1




MK-SOE1




MK-WOE1

H2




MK-NOE2




MK-SOE2


H3



MK-NOE3




MK-SOE3




MK-WOE3

H4




MK-NOE4




MK-SOE4


H5



MK-NOE5



MK-SOE5



MK-WOE5

OT-Equators System Components

OT-Equators are sold in kits



Impression coping cap



OE-IC002

Analog



OE-RSM60

Metal housing and
transfer caps



OE-MH001
OE-RCB01
OE-RCY01
OE-RCP01
OE-RCW01
OE-RCV01

Handling tool



ET-IT002

Driver for 0.05 inch hex.



MT-RDS30 - Short
MT-RDL30 - Long

Overdenture
Ball
Attachment System

Tightening torque: 30 Ncm

PLATFORM :


H1

SEVEN/M4
NP




MB-N2330

SEVEN/M4
SP



MB-N1375

SEVEN/M4
WP



MW-B1500

ANGULATED
SP


H2




MB-N3330



MB-N2375



MB-AN152




MB-AN252

H3



MB-N4330



MB-N3375



MW-B3500

H4



MB-N4375

H5



MW-B5500

Ball Attachment System Components

Analog




MB-RS375

Titanium ball cap



MB-TTP10

Plastic mounting sleeve



MB-PR010

Gold ball cap




MB-GPS10

Plastic ball cap




MB-SF200

Cap spring



MB-TBS10

Cap spring



MB-TBS20

Standard plastic cap for metal housing




MB-PPC10

Soft plastic cap for metal housing



MB-PPR10

Plastic disc



MB-DB235

Double metal housing



MB-DMH10

Ball attachment anchor key



MT-RB225

Prosthetic Screws and Drivers

Tightening torque: 20/30 Ncm

Prosthetic screw



MN-S0160 - NP
MD-S0220 - SP

Medium prosthetic
screw



MD-S0222

Long prosthetic screw



MD-S0224

Gold prosthetic screw



MD-G0220

Extra-short driver for 0.05
inch hex.



MT-RDM30

Short driver for 0.05 inch
hex.



MT-RDS30

Long driver for 0.05 inch
hex.




MT-RDL30

Extra-long driver for 0.05
inch hex.




MT-RDE30

Short motor adapter for
0.05 inch hex.



MT-SM005

Long motor adapter for
0.05 inch hex.



MT-LM005

Prosthetic Screws and Keys for Angled Screw Channel

Prosthetic screw for
angled screw channel*




MM-SA160

EZ-Base prosthetic screw



MN-E0160 - NP
MD-E0220 - SP




A unique screw and driver tip
design enable extreme angulation
and a wide insertion path.

Ratchet key for angled screw
channel



MT-ELR10 - Long
MT-ESR10 - Short


Motor key for angled screw
channel



MT-ELM10 - Long
MT-ESM10 - Short

MK-0039
Prosthetic Tools Kit





All rights reserved. No part of this publication may be reproduced, transcribed, stored in an electronic retrieval system, translated into any language or computer language, or be transmitted in any form whatsoever, without the prior written consent of the publisher. Warning: MIS products referred to in this document should be used by licensed dentists only.

Products may not be available in some regulated markets.
Please contact your local Dentsply Sirona sales office for current product assortment and availability.

